

Vernier Spectrophotometer User's Guide




1. Use a USB cable to connect a spectrometer *to the computer* (not to the USB hub)
2. Open Logger Pro
3. **Calibrate the spectrometer.**


- a. Prepare a *blank* by filling an empty cuvette $\frac{3}{4}$ full with DI water. Place the blank cuvette in the spectrometer with the *clear sides of the cuvette in the path of the light source*.
- b. Select **Calibrate Spectrometer** from the Experiment menu. The calibration dialog box will display the message: "Waiting seconds for lamp to warm up." The *minimum* warm up time is one minute, but best to wait ~3 minutes. Follow the instructions in the dialog box to complete the calibration. Click .

4. Absorption Spectrum and the Determination of the λ_{\max} of a colored solution

- a. Empty the blank cuvette and rinse it twice with small amounts of the colored solution. Fill the cuvette $\frac{3}{4}$ full with the solution and place it in the spectrometer.
- b. Click . A graph of the colored solution's absorption spectrum will be displayed. Note that one area of the graph contains a peak absorbance. Click to complete the analysis.
- c. *To save your graph of absorbance vs. wavelength:* Select **Store Latest Run** from the *Experiment menu*. If required by your instructor you may print the graph (absorption spectrum)—Print *two* copies of the graph—one for the white pages and one for the yellow pages of your lab report.





5. Configure the Spectrometer for Data Collection:

- a. Click the data collection icon, , on the toolbar. A dialog box will appear:
 - Select **Abs vs. Concentration** under *Set Collection Mode*. The wavelength of peak absorbance (λ_{\max}) will be automatically selected. If you wish to select a new wavelength, click on the graph or check the box next to the desired wavelength. Click to proceed.

6. Collection of absorbance-concentration data and the making a Calibration Curve (Beer's Law Plot)

- a. Place a cuvette containing the first standard solution in the spectrometer. Click . When the absorbance reading stabilizes, click (**Do not click Stop!!**). Enter the concentration of the standard solution and click (**Do not click Stop!!**). Record the absorbance of the solution.
- b. Discard the cuvette contents. Using the next standard solution, rinse and fill the cuvette $\frac{3}{4}$ full. Wipe the cuvette and place it in the spectrometer. When the absorbance reading stabilizes, click (**Do not click Stop!!**). Enter the concentration of Solution 3 (in moles/L) and click (**Do not click Stop!!**). Record the absorbance of solution 3 in your lab notebook.
- c. Repeat Step 5b for the remaining standard solutions.
- d. When you have finished testing the last of the standard solutions, click .



7. To determine the best-fit line equation for the standard solutions, click the linear fit button, , on the toolbar. *Write the equation of the standard curve in your lab notebook.*

- a. **Print** two copies of your absorbance vs. concentration graph, one copy for the white pages and the other for the yellow pages of your lab report.
- b. Select **Save As** from the File menu and save your experiment file.