

Barbie Bungee Jump

Lots of things vibrate or oscillate. A vibrating tuning fork, a moving child's playground swing, and the loudspeaker in a radio are all examples of physical vibrations. There are also electrical and acoustical vibrations, such as radio signals and the sound you get when blowing across the top of an open bottle.

One simple system that vibrates is a Barbie Doll hanging from a spring. The force applied by an ideal spring is proportional to how much it is stretched or compressed. Given this force behavior, the up and down motion of the Barbie Doll is called *simple harmonic* and the position can be modeled with

$$y = A\sin(2\pi ft + \phi)$$

In this equation, y is the vertical displacement from the equilibrium position, A is the amplitude of the motion, f is the frequency of the oscillation, t is the time, and ϕ is a phase constant. This experiment will clarify each of these terms.




Figure 1

OBJECTIVES

- Measure the position and velocity as a function of time for an oscillating Barbie Doll and spring system.
- Compare the observed motion of a Barbie Doll and spring system to a mathematical model of simple harmonic motion.
- Determine the amplitude, period, and phase constant of the observed simple harmonic motion.

MATERIALS






computer
Vernier computer interface
Logger Pro
Vernier Motion Detector
200 g and 300 g Barbie Dolls

ring stand, rod, and clamp
Rubber bands
twist ties
wire basket

PRELIMINARY QUESTIONS

1. Attach the Barbie Doll to a single rubber band and hold the free end of the rubber band in your hand, so the doll and spring hang down with the doll at rest. Lift the doll about 10 cm (4 inches) and release. Observe the motion. Sketch a graph of position vs. time for the Barbie Doll.
2. Just below the graph of position vs. time, and using the same length time scale, sketch a graph of velocity vs. time for the Barbie Doll.

PROCEDURE


1. Attach three rubber bands to a horizontal rod connected to the ring stand and hang the doll from the rubber bands as shown in Figure 1 using twist ties so the doll cannot fall.
2. Connect the Motion Detector to the DIG/SONIC 1 channel of the interface. If the Motion Detector has a switch, set it to Normal. 
3. Place the Motion Detector at least 75 cm below the doll. Make sure there are no objects near the path between the detector and Barbie Doll, such as a table edge.
4. Open the file “15 Simple Harmonic Motion” from the *Physics with Vernier* folder.
5. Make a preliminary run to make sure things are set up correctly. Lift the doll upward a few centimeters and release. The Barbie Doll should oscillate along a vertical line only. Click  to begin data collection.
6. After 10 s, data collection will stop. The position graph should show a clean sinusoidal curve. If it has flat regions or spikes, reposition the Motion Detector and try again.
7. Compare the position graph to your sketched prediction in the Preliminary Questions. How are the graphs similar? How are they different? Also, compare the velocity graph to your prediction.
8. Measure the equilibrium position of the Barbie Doll. Do this by allowing the Barbie Doll to hang free and at rest. Click  to begin data collection. After collection stops, click the Statistics button, , to determine the average distance from the detector. Record this distance (y_0) in your data table.
9. Now lift the Barbie Doll upward about 5 cm and release it. The Barbie Doll should oscillate along a vertical line only. Click  to collect data. Examine the graphs. The pattern you are observing is characteristic of simple harmonic motion.
10. Using the position graph, measure the time interval between maximum positions. This is the *period*, T , of the motion. The frequency, f , is the reciprocal of the period, $f = 1/T$. Based on your period measurement, calculate the frequency. Record the period and frequency of this motion in your data table.
11. The amplitude, A , of simple harmonic motion is the maximum distance from the equilibrium position. Estimate values for the amplitude from your position graph. Enter the values in your data table. If you drag the mouse from one peak to another you can read the Δx time interval.
12. Repeat Steps 8–11 with the same Barbie Doll, moving with a larger amplitude than in the first run. Record this data in run 2 of the data table.

13. Add 100 g mass to the Barbie Doll and repeat Steps 8–11. Use an amplitude of about 5 cm. Keep a good run made with this 100 g plus Barbie Doll on the screen. You will use it for several of the Analysis questions. Record this in Run 3 of the data table.

DATA TABLE

Run	Barbie (g)	y_0 (cm)	A (cm)	T (s)	f (Hz)
1					
2					
3					

ANALYSIS

- View the graphs of the last run on the screen. Compare the position vs. time and the velocity vs. time graphs. How are they the same? How are they different?
- Turn on the Examine mode by clicking the Examine button, . Move the mouse cursor back and forth across the graph to view the data values for the last run on the screen. Where is the Barbie Doll when the velocity is zero? Where is the Barbie Doll when the velocity is greatest?
- Does the frequency, f , appear to depend on the amplitude of the motion? Do you have enough data to draw a firm conclusion?
- Does the frequency, f , appear to depend on the Barbie Doll used? Did it change much in your tests?
- You can compare your experimental data to the sinusoidal function model using the Manual Curve Fit feature of *Logger Pro*. Try it with your 300 g data. The model equation in the introduction, which is similar to the one in many textbooks, gives the displacement from equilibrium. However, your Motion Detector reports the distance from the detector. To compare the model to your data, add the equilibrium distance to the model; that is, use

$$y = A \sin(2\pi ft + \phi) + y_0$$

where y_0 represents the equilibrium distance.

- Click once on the position graph to select it.
- Choose Curve Fit from the Analyze menu.
- Select Manual as the Fit Type.
- Select the Sine function from the General Equation list.
- The Sine equation is of the form $y = A \sin(Bt + C) + D$. Compare this to the form of the equation above to match variables; e.g., ϕ corresponds to C, and $2\pi f$ corresponds to B.
- Adjust the values for A, B and D to reflect your values for A, f and y_0 . You can either enter the values directly in the dialog box or you can use the up and down arrows to adjust the values.
- The phase parameter ϕ is called the *phase constant* and is used to adjust the y value reported by the model at $t = 0$ so that it matches your data. Since data collection did not necessarily begin when the Barbie Doll was at the equilibrium position, ϕ is needed to achieve a good match.

- h. The optimum value for ϕ will be between 0 and 2π . find a value for ϕ that makes the model come as close as possible to the data of your 300 g experiment. You may also want to adjust y_0 , A , and f to improve the fit. Write down the equation that best matches your data.
6. Does the model fit the data well? How can you tell?
7. Predict what would happen to the plot of the model if you doubled the parameter for A by sketching both the current model and the new model with doubled A . Now double the parameter for A in the manual fit dialog box to compare to your prediction.
8. Similarly, predict how the model plot would change if you doubled f , and then check by modifying the model definition.
9. Click , and optionally print your graph.

EXTENSIONS

1. Investigate how changing the spring amplitude changes the period of the motion. Take care not to use too large an amplitude so that the Barbie Doll does not come closer than 40 cm to the detector or fall from the spring.
2. How will *damping* change the data? Tape an index card to the top of the Barbie Doll and collect additional data. You may want to take data for more than 10 seconds. Does the model still fit well in this case?
3. Do additional experiments to discover the relationship between the Barbie Doll and the period of this motion.